

Hsp70 Rabbit mAb


ABZ5233

Product Overview


Name	Hsp70 Rabbit mAb
Accession(Primary)	P0DMV8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; FC-1:20; IP-1:20-1:50

Target information (P0DMV8)

Gene Name	HSPA1A
About Gene Symbol	



Western blot analysis of Hsp70 in HeLa, HCT116, A549 lysates using Hsp70 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Hsp70 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.