

C7 Rabbit mAb


ABZ5234

Product Overview


Name	C7 Rabbit mAb
Accession(Primary)	P10643
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 94 kDa; Observed MW: 94 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200

Target information (P10643)


Gene Name	C7
About Gene Symbol	



Western blot analysis of C7 in C6, PC-12, 3T3, Raw264.7 lysates using C7 antibody.



Western blot analysis of C7 in HeLa, 293 lysates using C7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using C7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.