

PSMD14 Rabbit mAb


ABZ5235

Product Overview


Name	PSMD14 Rabbit mAb
Accession(Primary)	O00487
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:1000

Target information (O00487)


Gene Name	PSMD14
About Gene Symbol	




Western blot analysis of PSMD14 in Raw264.7, 3T3 lysates using PSMD14 antibody.



Western blot analysis of PSMD14 in HeLa, COS7, PC-12 lysates using PSMD14 antibody.



Immunocytochemistry analysis of PSMD14(green)
in HeLa using PSMD14 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-
embedded Human Cholangiocarcinoma using
PSMD14 antibody. High-pressure and temperature
Sodium Citrate pH 6.0 was used for antigen
retrieval.