

NSE Rabbit mAb


ABZ5236

Product Overview


Name	NSE Rabbit mAb
Accession(Primary)	P09104
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:20-1:50; FC-1:20-1:100

Target information (P09104)


Gene Name	ENO2
About Gene Symbol	




Western blot analysis of NSE in mouse brain, rat brain lysates using NSE antibody.



Western blot analysis of Enolase2/NSE in Human Fetal Brain lysates using Enolase2/NSE antibody.



Immunocytochemistry analysis of NSE(green) in Hela using NSE antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using NSE antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.