

IKZF3 Rabbit mAb


ABZ5238

Product Overview


Name	IKZF3 Rabbit mAb
Accession(Primary)	Q9UKT9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100; FC-1:10-1:100

Target information (Q9UKT9)

Gene Name	IKZF3
About Gene Symbol	



Western blot analysis of IKZF3 in Raji, Ramos lysates using IKZF3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using IKZF3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.