

PARP1 Rabbit mAb


ABZ5240

Product Overview


Name	PARP1 Rabbit mAb
Accession(Primary)	P09874
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 113 kDa; Observed MW: 116,24 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; FC-1:20-1:50

Target information (P09874)


Gene Name	PARP1
About Gene Symbol	



Western blot analysis of PARP in 293, Jurkat lysates using PARP1 antibody.



Western blot analysis of PARP in HeLa lysates using PARP antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PARP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.