

# Lysozyme Rabbit mAb


ABZ5241

## Product Overview


|                             |                                                             |
|-----------------------------|-------------------------------------------------------------|
| <b>Name</b>                 | Lysozyme Rabbit mAb                                         |
| <b>Accession(Primary)</b>   | P61626                                                      |
| <b>Host Species</b>         | Rabbit                                                      |
| <b>Clonality</b>            | Monoclonal                                                  |
| <b>Isotype</b>              | IgG                                                         |
| <b>Species Reactivity</b>   | Human, Mouse                                                |
| <b>Calculated MW</b>        | Calculated MW: 17 kDa; Observed MW: 17 kDa                  |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:500-1:5000;<br>ICC/IF-1:100-1:200 |

## Target information (P61626)

|                          |     |
|--------------------------|-----|
| <b>Gene Name</b>         | LYZ |
| <b>About Gene Symbol</b> |     |



Western blot analysis of Lysozyme in HL-60 lysates using Lysozyme antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Lysozyme antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.