

Cytokeratin 8 Rabbit mAb


ABZ5242

Product Overview


Name	Cytokeratin 8 Rabbit mAb
Accession(Primary)	P05787
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:500; FC-1:20-1:50; IP-1:20

Target information (P05787)


Gene Name	KRT8
About Gene Symbol	



Western blot analysis of Keratin 8 in MCF-7, HeLa lysates using Cytokeratin 8 antibody.



Immunocytochemistry analysis of Cytokeratin 8(green) in HeLa using Cytokeratin 8 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Cytokeratin 8 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.