

Plastin L Rabbit mAb


ABZ5244

Product Overview


Name	Plastin L Rabbit mAb
Accession(Primary)	P13796
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; FC-1:200-1:1000; ICC/IF-1:50-1:100

Target information (P13796)

Gene Name	LCP1
About Gene Symbol	



Immunocytochemistry analysis of Plastin L (green) in Jurkat using Plastin L antibody, and DAPI(blue).



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Plastin L antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.