

CK18 Rabbit mAb


ABZ5245

Product Overview


Name	CK18 Rabbit mAb
Accession(Primary)	P05784
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:1000; ICC/IF-1:50-1:100

Target information (P05784)


Gene Name	Krt18
About Gene Symbol	



Western blot analysis of Keratin 18 in HeLa, K562 lysates using Cytokeratin 18 antibody.



Immunocytochemistry analysis of Cytokeratin 18(green) in HeLa using Cytokeratin 18 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Cytokeratin 18 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.