

VDAC1/Porin Rabbit mAb


ABZ5247

Product Overview


Name	VDAC1/Porin Rabbit mAb
Accession(Primary)	P21796
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:500-1:1000

Target information (P21796)


Gene Name	VDAC1
About Gene Symbol	



Western blot analysis of Calreticulin in (1) HepG2 lysates; (2) Jurkat lysates using VDAC1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using VDAC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of VDAC1 in C2C12, THP-1, 3T3, Jurakt and Hela lysates using VDAC1 antibody

