

Bcl 2 Rabbit mAb


ABZ5249

Product Overview


Name	Bcl 2 Rabbit mAb
Accession(Primary)	P10415
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; FC-1:50-1:100; ICC/IF-1:50-1:100

Target information (P10415)


Gene Name	BCL2
About Gene Symbol	BCL2 (BCL2 Apoptosis Regulator) is a protein phosphatase involved in apoptosis. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Bcl2 in 293T, 3T3 and HeLa lysates using Bcl2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Bcl2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Bcl2(green) in MCF-7 using Bcl2 antibody, and DAPI(blue)

