

Cyclin D1 Rabbit mAb


ABZ5250

Product Overview


Name	Cyclin D1 Rabbit mAb
Accession(Primary)	P24385
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:500; ICC/IF-1:20-1:50; IP-1:20-1:50

Target information (P24385)


Gene Name	CCND1
About Gene Symbol	CCND1 (Cyclin D1) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cyclin D1 in 3T3 lysates using Cyclin D1 antibody.



Western blot analysis of Cyclin D1 in MCF-7, C6, PC-12 lysates using Cyclin D1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Cyclin D1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.