

CD9 Rabbit mAb


ABZ5251

Product Overview


Name	CD9 Rabbit mAb
Accession(Primary)	P21926
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:100; FC-1:20; IP-1:20

Target information (P21926)


Gene Name	CD9
About Gene Symbol	



Western blot analysis of CD9 in U87-MG lysates using CD9 antibody.



Immunocytochemistry analysis of CD9(green) in Hela using CD9 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using CD9 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.