

## Fas Rabbit mAb


ABZ5253

### Product Overview

|                             |                                                            |
|-----------------------------|------------------------------------------------------------|
| <b>Name</b>                 | Fas Rabbit mAb                                             |
| <b>Accession(Primary)</b>   | P25445                                                     |
| <b>Host Species</b>         | Rabbit                                                     |
| <b>Clonality</b>            | Monoclonal                                                 |
| <b>Isotype</b>              | IgG                                                        |
| <b>Species Reactivity</b>   | Human                                                      |
| <b>Calculated MW</b>        | Calculated MW: 38 kDa; Observed MW: 48 kDa                 |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500;<br>ICC/IF-1:200-1:500 |

### Target information (P25445)

|                          |     |
|--------------------------|-----|
| <b>Gene Name</b>         | FAS |
| <b>About Gene Symbol</b> |     |



Western blot analysis of Fas in Ramos lysates  
using Fas antibody.