

# Hydroxyl-Histone H2A (Tyr39) Rabbit mAb

ABZ5258


## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Hydroxyl-Histone H2A (Tyr39) Rabbit mAb    |
| <b>Accession(Primary)</b>   | P04908                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 14 kDa; Observed MW: 14 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:500-1:2000          |


## Target information (P04908)

**Gene Name** H2AC4

**About Gene Symbol**



Western blot analysis of Calreticulin in (1) NIH/3T3 lysates; (2) A549 lysates using Hydroxyl-Histone H2A (Tyr39) antibody.



Immunohistochemistry analysis of paraffin-embedded mouse colon using Histone H2A (Hydroxyl- Y39) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.