

Phospho-Syk (Tyr323) Rabbit mAb


ABZ5263

Product Overview

Name	Phospho-Syk (Tyr323) Rabbit mAb
Accession(Primary)	P43405
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 72 kDa; Observed MW: 72 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:100-1:200; ICC/IF-1:100-1:200; FC-1:20

Target information (P43405)

Gene Name	SYK
About Gene Symbol	



Western blot analysis of Phospho-Syk (Tyr323) in K562 lysates using Phospho-Syk (Tyr323) antibody.