

## Phospho-STAT1 (Ser727) Rabbit mAb


ABZ5266

### Product Overview


|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Phospho-STAT1 (Ser727) Rabbit mAb          |
| <b>Accession(Primary)</b>   | P42224                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 87 kDa; Observed MW: 87 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200        |

### Target information (P42224)

|                          |                                                                                                                                                                                                                                                                    |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | STAT1                                                                                                                                                                                                                                                              |
| <b>About Gene Symbol</b> | STAT1 (Signal Transducer And Activator Of Transcription 1)<br>is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cell |



Western blot analysis of Phospho-STAT1 (S727) in HeLa lysates treated with etoposide using Phospho-STAT1 (Ser727) antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Phospho-STAT1 (S727) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.