


Phospho-IRE1 (Ser724) Rabbit mAb

ABZ5270

Product Overview

Name	Phospho-IRE1 (Ser724) Rabbit mAb
Accession(Primary)	O75460
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 110 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (O75460)

Gene Name ERN1**About Gene Symbol**

Western blot analysis of Phospho-IRE1 (S724) in (1) K562 lysates treated with AP; (2) Untreated K562 lysates using Phospho-IRE1 (Ser724) antibody.

