

Phospho-LRRK2 (Ser395) Rabbit mAb

ABZ5273


Product Overview

Name	Phospho-LRRK2 (Ser395) Rabbit mAb
Accession(Primary)	Q5S007
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 286 kDa; Observed MW: 286 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q5S007)

Gene Name LRRK2

About Gene Symbol



Western blot analysis of Filamin A
Phosphorylation in WTLRRK2 lysates treated LRRK2
using Phospho-LRRK2 (Ser395) antibody.

