

Phospho-HSF1 (Ser326) Rabbit mAb


ABZ5274

Product Overview

Name	Phospho-HSF1 (Ser326) Rabbit mAb
Accession(Primary)	Q00613
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 57 kDa; Observed MW: 82 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q00613)

Gene Name	HSF1
About Gene Symbol	HSF1 (Heat Shock Transcription Factor 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-HSF1 (S326) in HeLa lysates treated with heat using Phospho-HSF1 (Ser326) antibody.