

Phospho-c Myb (Ser11) Rabbit mAb


ABZ5277

Product Overview

Name	Phospho-c Myb (Ser11) Rabbit mAb
Accession(Primary)	P10242
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 72 kDa; Observed MW: 72,80 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (P10242)

Gene Name	MYB
About Gene Symbol	



Western blot analysis of Phospho-vMyb + c-Myb (S11) in Ramos lysates using Phospho-c Myb (Ser11) antibody.