

Phospho-Histone H2A (Ser129) Rabbit mAb

ABZ5281


Product Overview

Name	Phospho-Histone H2A (Ser129) Rabbit mAb
Accession(Primary)	P04912
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P04912)

Gene Name HTA2

About Gene Symbol



Western blot analysis of Phospho-Histone H2A (S129) in *Saccharomyces cerevisiae* lysates treated with Methyl methanesulfonate using Phospho-Histone H2A (Ser129) antibody.

