

Phospho-Synapsin I (Ser9) Rabbit mAb

ABZ5284


Product Overview

Name	Phospho-Synapsin I (Ser9) Rabbit mAb
Accession(Primary)	P17600
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 74 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P17600)

Gene Name SYN1

About Gene Symbol



Western blot analysis of Phospho-Synapsin I (S9) in (1) Human brain lysates; (2) Human brain lysates treated with AP using Phospho-Synapsin I (Ser9) antibody.

