

Phospho-BRAF (Thr401) Rabbit mAb

ABZ5289


Product Overview

Name	Phospho-BRAF (Thr401) Rabbit mAb
Accession(Primary)	P15056
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P15056)

Gene Name BRAF

About Gene Symbol



Western blot analysis of Phospho-B Raf (T401) in (1) PC-12 lysates; (2) PC-12 lysates treated with TPA using Phospho-BRAF (Thr401) antibody.

