

Phospho-c-Jun (Ser63) Rabbit mAb


ABZ5290

Product Overview


Name	Phospho-c-Jun (Ser63) Rabbit mAb
Accession(Primary)	P05412
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (P05412)


Gene Name	JUN
About Gene Symbol	JUN (Jun Proto-Oncogene, AP-1 Transcription Factor Subunit) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of c-Jun Phosphorylation in NIH/3T3 lysates treated with Anisomycin using Phospho-c-Jun (Ser63) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma using Phospho-c-Jun (S63) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Phospho-c-Jun (Ser63) in HeLa treated with anisomycin using Phospho-c-Jun (S63) antibody.