

Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) Rabbit mAb


ABZ5291

Product Overview

Name	Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) Rabbit mAb
Accession(Primary)	O96013
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 64,80,75 kDa; Observed MW: 64,80,75 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:50-1:100

Target information (O96013)

Gene Name	PAK4/5/6
About Gene Symbol	



Western blot analysis of Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602) in U87-MG lysates using Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) antibody.