

Phospho-SIRT1 (Ser47) Rabbit mAb

ABZ5293


Product Overview

Name	Phospho-SIRT1 (Ser47) Rabbit mAb
Accession(Primary)	Q96EB6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 82 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (Q96EB6)

Gene Name SIRT1

About Gene Symbol



Western blot analysis of Phospho-SIRT1(Ser47) in A549 lysates using Phospho-SIRT1 (Ser47) antibody.

