

Phospho-160 kD Neurofilament Medium (Ser614/Ser619)


ABZ5294

Product Overview

Name	Phospho-160 kD Neurofilament Medium (Ser614/Ser619) Rabbit mAb
Accession(Primary)	P07197
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 102 kDa; Observed MW: 160 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (P07197)

Gene Name	NEFM
About Gene Symbol	NEFM (Neurofilament Medium Chain) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-Neurofilament Medium (Ser614/619) in K562, Ramos, Daudi lysates using Phospho-Neurofilament Medium (Ser614/Ser619) antibody.