

Phospho-PRKAR2A (Ser99) Rabbit mAb


ABZ5296

Product Overview

Name	Phospho-PRKAR2A (Ser99) Rabbit mAb
Accession(Primary)	P13861
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 51 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IF-1:200-1:500

Target information (P13861)

Gene Name	PRKAR2A
About Gene Symbol	PRKAR2A (Protein Kinase CAMP-Dependent Type II Regulatory Subunit Alpha) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research



Western blot analysis of Phospho-PKA R2 (Ser99)
in K562 lysates using Phospho-PKA RII alpha
(Ser99) antibody.