

MonoMethyl-Histone H2B (Arg79) Rabbit mAb

ABZ5302


Product Overview

Name	MonoMethyl-Histone H2B (Arg79) Rabbit mAb
Accession(Primary)	Q16778
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (Q16778)

Gene Name H2BC21

About Gene Symbol



Western blot analysis of Histone H2B (mono methyl R79) in HeLa lysates using MonoMethyl-Histone H2B (Arg79) antibody.

