

[KD]IRF3 Rabbit mAb


ABZ5306

Product Overview


Name	[KD]IRF3 Rabbit mAb
Accession(Primary)	Q14653
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 47 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100; FC-1:100-1:200

Target information (Q14653)


Gene Name	IRF3
About Gene Symbol	IRF3 (Interferon Regulatory Factor 3) is an interferon regulatory transcription factor involved in JAK-STAT/cytokine signaling and innate immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of IRF3 in (1) HeLa lysates;
(2) 3T3 lysates using IRF3 antibody.



Immunohistochemistry analysis of paraffin-
embedded Human cervix carcinoma using IRF3
antibody. High-pressure and temperature Sodium
Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of IRF3 in Jurkat using IRF3 antibody.