

[KD]Cyclophilin F Rabbit mAb


ABZ5311

Product Overview

Name	[KD]Cyclophilin F Rabbit mAb
Accession(Primary)	P30405
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:5000; ICC/IF-1:50-1:100; FC-1:200-1:500

Target information (P30405)

Gene Name	PPIF
About Gene Symbol	PPIF (Peptidylprolyl Isomerase F) is a peptidylprolyl isomerase (cyclophilin) involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cyclophilin F in HeLa
lysates using Cyclophilin F antibody.