

[KO]Rab9 Rabbit mAb


ABZ5315

Product Overview

Name	[KO]Rab9 Rabbit mAb
Accession(Primary)	P51151
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200

Target information (P51151)

Gene Name	RAB9A
About Gene Symbol	RAB9A (RAB9A, Member RAS Oncogene Family) is a small GTPase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Rab9 in HepG2 lysates
using Rab9 antibody.