

[KO]PARK7/DJ1 Rabbit mAb


ABZ5316

Product Overview


Name	[KO]PARK7/DJ1 Rabbit mAb
Accession(Primary)	Q99497
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 20 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:500; FC-1:20-1:50; IP-1:20

Target information (Q99497)

Gene Name	PARK7
About Gene Symbol	PARK7 (Parkinsonism Associated Deglycase) is a metabolic enzyme involved in cellular metabolism. It is an important research target in cancer research and neurodegenerative disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PARK7 in HeLa lysates using PARK7/DJ1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human bladder cancer using PARK7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.