

[KO]Histone H1 Rabbit mAb


ABZ5317

Product Overview


Name	[KO]Histone H1 Rabbit mAb
Accession(Primary)	P07305
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (P07305)


Gene Name	H1-0
About Gene Symbol	H1-0 (H1.0 Linker Histone) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Histone H1.0 in Human kidney lysates using Histone H1.0 antibody.



Immunohistochemistry analysis of paraffin-embedded Human testis using Histone H1.0 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Histone H1.0 in HepG2 using Histone H1.0 antibody.

