

[KD]A-Raf Rabbit mAb


ABZ5318

Product Overview

Name	[KD]A-Raf Rabbit mAb
Accession(Primary)	P10398
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (P10398)

Gene Name	ARAF
About Gene Symbol	ARAF (A-Raf Proto-Oncogene, Serine/Threonine Kinase) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ARAF in A375 lysates
using A RAF antibody.