

[KD]4E BP1 Rabbit mAb


ABZ5319

Product Overview


Name	[KD]4E BP1 Rabbit mAb
Accession(Primary)	Q13541
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 13 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:20-1:50; IP-1:20

Target information (Q13541)


Gene Name	EIF4EBP1
About Gene Symbol	EIF4EBP1 (Eukaryotic Translation Initiation Factor 4E Binding Protein 1) is a protein kinase and signaling regulator involved in mTOR/PI3K/Akt signaling and translational regulation. It is an important research target in cancer, cell signaling, and therap



Western blot analysis of eIF4EBP1 in K562 lysates using 4E BP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma using eIF4EBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of 4E BP1 in HeLa using eIF4EBP1 antibody.

