

[KD]SDHA Rabbit mAb


ABZ5320

Product Overview

Name	[KD]SDHA Rabbit mAb
Accession(Primary)	P31040
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:500; ICC/IF-1:100-1:200; FC-1:20; IP-1:10-1:50

Target information (P31040)

Gene Name	SDHA
About Gene Symbol	



Western blot analysis of SDHA in (1) HeLa lysates; (2) mouse brain lysates; (3) rat brain lysates; using SDHA antibody.