


[KD]Coilin Rabbit mAb

ABZ5325

Product Overview

Name	[KD]Coilin Rabbit mAb
Accession(Primary)	P38432
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 63 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:500

Target information (P38432)

Gene Name COIL**About Gene Symbol**

Western blot analysis of Coilin in Jurkat lysates using Coilin antibody.

