

[KD]Integrin alpha V Rabbit mAb


ABZ5326

Product Overview

Name	[KD]Integrin alpha V Rabbit mAb
Accession(Primary)	P06756
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 116 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; IP-1:20-1:50

Target information (P06756)

Gene Name	ITGAV
About Gene Symbol	ITGAV (Integrin Subunit Alpha V) is an integrin cell adhesion receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Integrin alpha V in (1) A549 lysates; (2) mouse spleen lysates; (3) rat brain lysates using Integrin alpha V antibody.