

[KD]Aurora B Rabbit mAb


ABZ5327

Product Overview


Name	[KD]Aurora B Rabbit mAb
Accession(Primary)	Q96GD4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:20; FC-1:200-1:500; IP-1:20

Target information (Q96GD4)


Gene Name	AURKB
About Gene Symbol	



Western blot analysis of Aurora B in HeLa lysates using Aurora B antibody.



Western blot analysis of Aurora B in Jurkat and HeLa lysates using Aurora B antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil, using Aurora B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.