

[KO]CD166 Rabbit mAb


ABZ5329

Product Overview

Name	[KO]CD166 Rabbit mAb
Accession(Primary)	Q13740
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50; ICC/IF-1:50-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (Q13740)

Gene Name	ALCAM
About Gene Symbol	ALCAM (Activated Leukocyte Cell Adhesion Molecule) is a cell surface marker involved in adaptive immunity. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CD166 in Daudi lysates
using CD166 antibody.