

[KO]FKBP51 Rabbit mAb


ABZ5330

Product Overview

Name	[KO]FKBP51 Rabbit mAb
Accession(Primary)	Q13451
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 51 kDa; Observed MW: 51 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20; FC-1:50-1:100

Target information (Q13451)

Gene Name	FKBP5
About Gene Symbol	



Western blot analysis of FKBP51 in Caco2 lysates
using FKBP51 antibody.