

# [KO]Aryl Hydrocarbon Receptor Rabbit mAb


ABZ5335

## Product Overview

|                             |                                             |
|-----------------------------|---------------------------------------------|
| <b>Name</b>                 | [KO]Aryl Hydrocarbon Receptor Rabbit mAb    |
| <b>Accession(Primary)</b>   | P35869                                      |
| <b>Host Species</b>         | Rabbit                                      |
| <b>Clonality</b>            | Monoclonal                                  |
| <b>Isotype</b>              | IgG                                         |
| <b>Species Reactivity</b>   | Human                                       |
| <b>Calculated MW</b>        | Calculated MW: 96 kDa; Observed MW: 105 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; ICC/IF-1:100-1:200         |

## Target information (P35869)

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | AHR                                                                                                                                                                                                                                                                                |
| <b>About Gene Symbol</b> | AHR (Aryl Hydrocarbon Receptor) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of AhR in HEK293 lysates  
using Aryl Hydrocarbon Receptor antibody.