

[KD]CTNNA1 Rabbit mAb


ABZ5336

Product Overview


Name	[KD]CTNNA1 Rabbit mAb
Accession(Primary)	P35221
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 100 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:100; IP-1:50-1:100

Target information (P35221)

Gene Name	CTNNA1
About Gene Symbol	CTNNA1 (Catenin Alpha 1) is a cell junction/adhesion protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Catenin alpha 1 in HeLa lysates using CTNNA1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Catenin alpha 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.