

## [KD]DDX5 Rabbit mAb


ABZ5337

### Product Overview

|                             |                                                                                           |
|-----------------------------|-------------------------------------------------------------------------------------------|
| <b>Name</b>                 | [KD]DDX5 Rabbit mAb                                                                       |
| <b>Accession(Primary)</b>   | P17844                                                                                    |
| <b>Host Species</b>         | Rabbit                                                                                    |
| <b>Clonality</b>            | Monoclonal                                                                                |
| <b>Isotype</b>              | IgG                                                                                       |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                                         |
| <b>Calculated MW</b>        | Calculated MW: 69 kDa; Observed MW: 69 kDa                                                |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500;<br>ICC/IF-1:100-1:200; FC-1:100-1:200; IP-1:10-1:100 |

### Target information (P17844)

|                          |                                                                                                                                                                                                                                                                           |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | DDX5                                                                                                                                                                                                                                                                      |
| <b>About Gene Symbol</b> | DDX5 (DEAD-Box Helicase 5) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of DDX5 in HeLa lysates  
using DDX5 antibody.