

[KD]SERCA2 ATPase Rabbit mAb


ABZ5338

Product Overview

Name	[KD]SERCA2 ATPase Rabbit mAb
Accession(Primary)	P16615
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 115 kDa; Observed MW: 110,140 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:100-1:1000

Target information (P16615)

Gene Name	ATP2A2
About Gene Symbol	ATP2A2 (ATPase Sarcoplasmic/Endoplasmic Reticulum Ca ²⁺ Transporting 2) is a membrane transporter involved in membrane transport. It is an important research target in cardiovascular disease and is widely studied to elucidate cellular mechanisms and identi



Western blot analysis of SERCA2 in HeLa lysates
using SERCA2 ATPase antibody.