

[KD]Menin Rabbit mAb


ABZ5343

Product Overview

Name	[KD]Menin Rabbit mAb
Accession(Primary)	O00255
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:500; ICC/IF-1:20-1:50; IP-1:10-1:100

Target information (O00255)

Gene Name	MEN1
About Gene Symbol	MEN1 (Menin 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Menin in 293T lysates
using Menin antibody.