

ACTA1 (2B10) Mouse mAb

ABZ3926

Product Overview

Name	ACTA1 (2B10) Mouse mAb
Accession(Primary)	P68133
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (P68133)

Gene Name	ACTA1
About Gene Symbol	

